



Sterile Alcohol Prep Pads

For Optimal Antibacterial Action

- For topical cleansing prior to injection or venipuncture and manufactured according to current FDA regulations, GMP and USP requirements
- Provides rapid-acting skin antiseptics in a sterile formula
- Sterile alcohol prep pads are saturated with 70% isopropyl alcohol
- Soft and absorbent, non-woven pad
- Gamma sterilized and not made with natural rubber latex



Reorder No.	Size	Length Folded	Length Unfolded	Width	Quantity
1113	Medium	1.18"	2.36"	1.1"	10/200/Cs
1114	Medium	1.18"	2.69"	1.1"	20/100/Cs
1116	Large	1.77"	3.54"	1.18"	10/100/Cs

*Sizes will vary slightly based on product specification tolerances
 HCPCS: A4245

Find us on:



For more information contact your local Omni Health Solution sales representative OR Call us at: (888) 914-OMNI
info@omni-healthsolution.com