


Sienna™

COVID-19

ANTIGEN RAPID TEST CASSETTE (NASOPHARYNGEAL SWAB)



This rapid chromatographic immunoassay was created for the qualitative detection of COVID-19 antigen. The identification is based on the monoclonal antibodies specific for the nucleocapsid protein of SARS-CoV-2. It is intended to aid in the rapid diagnosis of COVID-19 infections.

The test is intended for professional in vitro diagnostic use. The kit includes 1 IFU, 25 individually pouched cassettes, 25 buffer tubes and 25 sterile swabs. Store as packaged at room temperature or refrigerated 2-30 degrees C / 36-86 degrees F. **DO NOT FREEZE.**




Dependable, accurate results

- Sensitivity: 93.8%
- Specificity: 99.2%
- LOD: 4.0×10^3 TCID₅₀/ml

Unmatched convenience

- Individual buffer vials: helping multiple operators at the same time
- Visual results: simple to interpret
- 10 minute test time: timely results



Omni Health Solution, CA LLC mission is bringing a range of product solutions to improve the health and wellness of people around the world.

The company's combined experience as C-suite executives at the world's largest health systems, experts in the high complexity lab space, healthcare entrepreneurs, and healthcare legal experts brings the expertise required to identify good science and solid business practices.
www.info@omni-healthsolution.com